


UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM.

**Current Status of the Relay Ascending Aortic PS-IDE:
What makes it different?**



Adam W. Beck, MD
Veith Symposium
November 2024
NY, New York

UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM DIVISION OF VASCULAR SURGERY & ENDOVASCULAR THERAPY © UAB. All Rights Reserved.

Disclosures

All proceeds to UAB

Consultant for:
Artivion, Cook Medical, Philips,
Medtronic, Terumo Aortic

Research support: Cook Medical, Endospin,
Medtronic, Philips, Terumo, W.L. Gore
& Associates

UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM DIVISION OF VASCULAR SURGERY & ENDOVASCULAR THERAPY © UAB. All Rights Reserved.

Additional Disclosures

- A PS-IDE using the Modified RelayPro Ascending Device is being reviewed by FDA through a pre-submission
- The Modified RelayPro Ascending Device is currently not available in the United States.
- An Ascending Relay Device is available in Europe through the custom-made program
- (Custom made devices are specifically made in accordance with a written prescription of any person authorized by national law by virtue of that person's professional qualifications; which gives (1) specific design characteristics provided under that person's responsibility and (2) is intended for the sole use of a particular patient exclusively to meet their individual conditions and needs. Custom made devices are not available in the US and availability is subject to local regulatory approval.)

UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM DIVISION OF VASCULAR SURGERY & ENDOVASCULAR THERAPY © UAB. All Rights Reserved.

Relay® Pro Platform

Low Profile

- 21FR-23FR OD

Appropriate Usable Length

- Long outer sheath (90cm working length)
- Ability to reach Zone 0



UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM DIVISION OF VASCULAR SURGERY & ENDOVASCULAR THERAPY © UAB. All Rights Reserved.

Relay® Pro Platform

Dual Sheath Technology

- Nitinol pre-curved inner catheter
- Inner sheath
- Outer sheath

Core Technologies

- Flared constraining sleeve
- Proximal clamping (allows re-positioning)
- Deployment sequence

Mechanical Advantage

- Controlled stent-graft expansion



UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM DIVISION OF VASCULAR SURGERY & ENDOVASCULAR THERAPY © UAB. All Rights Reserved.

Relay® Pro: Flexible Inner Sheath

Atraumatic navigation through the aortic arch

- Reduces potential trauma
- Staged deployment

Pre-curved inner catheter facilitates self-alignment to primary curve of the aorta


- Conforms to the ascending aorta

Enhanced constraining sleeve

- Optimized weave pattern allows for a reduction in profile
- Expands 9mm when exiting the outer sheath

Inner Sheath Flared End (V-patch)

- Allows for stent-graft expansion towards the inner curve for a correct apposition on the inner curve



UAB THE UNIVERSITY OF ALABAMA AT BIRMINGHAM DIVISION OF VASCULAR SURGERY & ENDOVASCULAR THERAPY © UAB. All Rights Reserved.

Modified Relay® Pro Platform

- Landing zone requirements:** The stent-graft must land 10mm proximal to the entry tear without compromising the function of the aortic root anatomy
- Non-bare stent proximally & distally
- 65mm & 85mm length
- Diameters from 32mm – 46mm
- No S-bar**

Stent-Graft Diameter (mm)	Infrarenal Proximal Vessel Diameter (mm)	Delivery System Fibresh-Box (D.B.)	
		65mm Length	85mm Length
32	28-32	21	21
34	30-31	21	21
36	30-31	22	22
38	34	22	22
40	35-36	22	22
42	37-38	23	23
44	39-40	23	23
46	41-42	23	23

Non-Bare Stent Configuration

65mm Length

85mm Length

UAB UNIVERSITY OF ALABAMA AT BIRMINGHAM


© UAB All Rights Reserved.

Study Description

- Aortic Team:** (Vascular + Cardiac Surgery + Cardiac Anesthesia)
- High and Prohibitive Risk for Open Surgery
- Zone 0 – Zone 0 deployment zone
- Indications:**
 - Acute and Chronic Ascending Aortic Dissection
 - Pseudoaneurysm
 - PAU/IMH

UAB UNIVERSITY OF ALABAMA AT BIRMINGHAM

© UAB All Rights Reserved.



Thank You

- awbeck@uabmc.edu
- <https://www.uab.edu/medicine/surgery/vascular-endovascular>

UAB UNIVERSITY OF ALABAMA AT BIRMINGHAM

DEPARTMENT OF SURGERY - DIVISION OF VASCULAR AND ENDOVASCULAR SURGERY

© UAB All Rights Reserved.