

Endovascular Thermal Septotomy Of The Dissection Flap With An Electrocautery Activated Guidewire To Treat Chronic TBADs:

When And How Does It Work: Complications And When It Should Be Avoided: Technical Tips And How To Avoid Pitfalls

Carlos H. Timaran, MD
 Professor of Surgery
 Chief, Endovascular Surgery
 Sam H. Phillips Jr MD Distinguished Chair in Surgery
 University of Texas Southwestern Medical Center - Dallas, TX

Division of Vascular and Endovascular Surgery

UTSouthwestern Medical Center

Transcatheter Electrosurgical Septotomy

A Disruptive Technique

Carlos H. Timaran, MD
 Professor of Surgery
 Chief, Endovascular Surgery
 Sam H. Phillips Jr MD Distinguished Chair in Surgery
 University of Texas Southwestern Medical Center - Dallas, TX

Division of Vascular and Endovascular Surgery

UTSouthwestern Medical Center

Disclosures

- Cook Medical Inc.
 - Consulting
 - Research support
 - Planning & proctoring
- Some devices presented here are investigational and have not been approved by the FDA
- Acknowledgement
 - Gustavo Oderich, MD
 Professor and Chief of Vascular Surgery
 UT Health & Memorial Hermann Texas Medical Center ,
 Houston, TX, USA

UTSouthwestern Medical Center

False lumen management




Department of Surgery

UTSouthwestern Medical Center

Transcatheter electrosurgery

2015

Structural heart disease:
BASILICA and LAMPOON
 techniques
(Khan et al. JACC 2018)




Division of Vascular and Endovascular Surgery

UTSouthwestern Medical Center

Transcatheter electrosurgery

2015

Structural heart disease:
BASILICA and LAMPOON
 techniques
(Khan et al. JACC 2018)



2019

Chronic dissection:
Transcatheter electrosurgical fenestration
(Kabbani L et al. JVSC 2023)



Transcatheter electrosurgical aortic septotomy (TEAS)

Division of Vascular and Endovascular Surgery

UTSouthwestern Medical Center

THE PRESENT AND FUTURE
JACC STATE-OF-THE-ART REVIEW
Transcatheter Electrosurgery
JACC State-of-the-Art Review

Jaffer M. Khan, BM BCh, PhD^{1,2}; Toby Rogers, BM BCh, PhD^{3,4}; Adam B. Cowanham, MD⁵; Vasilis C. Babalero, MD⁶; Dennis Korel Yikraz, MS⁷; Christopher G. Broce, MB ChB⁸; Daniel A. Herzig, PhD⁹; William H. Schenke, BS¹⁰; Kazuhiko Katsuyuki, MD¹¹; Robert J. Leshem, MD¹²

JACC:2020;75:1455-70

Department of Surgery **UT Southwestern Medical Center**

Transcatheter Electrosurgical Septotomy

ASTATO XS 20 300cm

High tension Stainless Steel Core
 Hydraulic Coating USP Class II (ASTATO)
 PTFE Coating

NAVICROSS® 0.018" Catheter
 Minimum Sheath
 Compatibility: 2.6 Fr

Denuding PTFE coating of Astato wire

UT Southwestern Medical Center

UT Southwestern Medical Center

Transcatheter Electrosurgical Aortic Septotomy



Division of Vascular and Endovascular Surgery **UT Southwestern Medical Center**

Transcatheter Electrosurgical Aortic Septotomy

Division of Vascular and Endovascular Surgery **UT Southwestern Medical Center**

Setup Electrocautery




Division of Vascular and Endovascular Surgery **UT Southwestern Medical Center**



Case report


- 70yo male presented with a thoracic aortic aneurysm
- First noted in 2012, managed conservatively
- March 2024, CTA revealed an increase in aneurysm sac to 61 mm and Type B dissection

Division of Vascular and Endovascular Surgery UT Southwestern Medical Center



Treatment Plan


- TEVAR using Thoracic Branch Endoprosthesis (TBE)
- Transcatheter Electrosurgical Aortic Septotomy




Division of Vascular and Endovascular Surgery | UTSouthwestern Medical Center

3-months later

- Presented to OSH with syncope, chest, and abdominal pain radiating to the back
- The patient was intubated and transferred to our institution




Division of Vascular and Endovascular Surgery | UTSouthwestern Medical Center








Plan

- Urgent endovascular thoracic aneurysm repair with TBE & left subclavian stenting
- Transcatheter Electrosurgical Aortic Septotomy




Division of Vascular and Endovascular Surgery | UTSouthwestern Medical Center





Hospital course

- Left hemothorax evacuation on POD 3
- Chest tubes removed on POD 5
- Discharged home on POD 8



Division of Vascular and Endovascular Surgery



Indications, safety, and effectiveness of transcatheter electrocatheter septotomy during endovascular repair of aortic dissections

Lucas Ruitter Kanamori, MD,* Emanuel R. Tenorio, MD, PhD,* Dora Babocs, MD,* Safa Savadi, MD,* Aidin Baghban-Oskoui, MD,* Ying Huang, MD, PhD,* Andres Figueroa, MD,* Mira Tamerbaum, MD,* Jose Eduardo Costa Filho, MD,* Mirza Baig, MD,* Thais A. Macedo, MD,* Carlos H. Timman, MD,* and Gustavo S. Oderich, MD,* Houston and Dallas, TX J Vasc Surg 2024;80:1396-406

UTHealth Houston McGovern Medical School UTSouthwestern Medical Center

2-Center Experience

	n = 36	Percent or IQR
Acute complicated dissection	3	8
Chronic post-dissection aneurysm	32	89
Aortic arch (Zone 0-3)	7	19
Thoracic aorta (Zone 4-5)	7	19
Thoracoabdominal aorta (Zone 4-9)	19	53
Infrarenal aorta and iliacs (Zone 9-11)	3	8
Maximum aortic diameter (mm)	60	52 - 70
Prior aortic repair	27	75
Prior open surgical repair	23	64
Prior endovascular aortic repair	15	42
Symptomatic/ ruptured aneurysm	10	28
Family History of aortic disease	3	8





#UTHealth Houston UTSouthwestern Medical Center

Extent of repair and procedural data

	n = 36	Percent or Mean
Arch branch stent graft	7	19
TEVAR (±Petticoat)	8	22
TAAA FB-EVAR	18	50
EVAR	3	8
Procedural data	11	31
Total operating time (min, mean ± SD)		335 ± 22
Total endovascular time (min, mean ± SD)		248 ± 162
Total fluoroscopy time (min, mean ± SD)		102 ± 67
Total contrast volume (ml, mean ± SD)		177 ± 56
Total Cumulative Air Kerma (Gy, mean ± SD)	36	2.0 ± 1.5

TEAS indication, extent and technical success

	n = 36	%
Indications		
True lumen compression (≤ 16mm)	28	78
Target artery from different lumen	19	53
Proximal or distal landing zone optimization	12	33
Organ or limb malperfusion	4	11
Technical success		
All patients (n = 36)	33	92
Chronic dissections (n = 33)	32	97
Acute dissections (n = 3)	1	33
Reasons for technical failure		
Dislodgement of dissection lamella, 2 (acute dissection)		
Inadvertent SMA dissection, 1 (chronic dissection)		
No arterial disruption or rupture		

Always look



Division of Vascular and Endovascular Surgery

UT Southwestern Medical Center

Conclusion

- TEAS can be performed safely and achieve a successful distal seal with false lumen exclusion in patients with **chronic** dissections
- TEAS facilitates endovascular repair, particularly in narrow true lumen and when target vessels arise from false lumen
- Full fluoroscopic visualization is needed to avoid injury
- Mural thrombus or thick septum should be avoided
- Longer follow up is needed to assess durability and frequency of endoleaks after FBEVAR for PD-TAAAs



Division of Vascular and Endovascular Surgery

UT Southwestern Medical Center