

**Remote Measurement Of Blood Flow In Any Artery Or Graft Using An Implantable Ultrasonic Sensor With Transmission To The Internet (From IntelVasc): How Does It Work, And Experimental Data**

**JEAN BISMUTH, MD**  
 Chief Division of Vascular Surgery, USF Health Morsani School of Medicine  
 Professor of Surgery University of South Florida  
 Adjunct Professor Rice University Institute of Biosciences and Bioengineering  
 Adjunct Professor Rice University Institute of Mechanical Engineering

TGH Tampa General Hospital | USF Health Morsani College of Medicine | TGH HEART & VASCULAR INSTITUTE

## Disclosures

WL Gore - speaker  
 Cook Medical - speaker  
 IntelVasc - CMO and co-founder

TGH HEART & VASCULAR INSTITUTE | USF Health

2 BEST PRACTICES



**INTELVASC**™  
 VASCULAR INSIGHTS

MONITORING VASCULAR OCCLUSION REMOTELY

### THE PROBLEM

**PERIPHERAL ARTERIAL DISEASE to USD...**

33M suffer from PAD in US | \$20B healthcare expenditure

**TIME TO DIAGNOSIS / COST**

EARLY	\$7,000
MODERATE	\$15,000
LATE	\$45,000

1.3M Surgeries a year


**Time Matters!**

30% PAD patients will have STROKE, HEART ATTACK, AMPUTATION

90% PAD patients are **ASYMPTOMATIC**

65% PAD patients miss MEDICAL CHECK-UP

### IDENTIFYING CLINICAL COMPLICATIONS EARLY




- MEDICAL VISIT**  
Patient attends an outpatient medical visit
- BIOSENSOR INSERTION**  
The specialist inserts the sensor subcutaneously
- AT-HOME MONITORING**  
Biosensor sends data to mobile app comfortably from home
- MOBILE BASED CLINICAL UPLOAD**  
Assessment data is uploaded to the cloud
- REMOTE CLINICAL EVALUATION**  
If vascular problems are detected, clinical visit is scheduled

**HIGHER PATIENT COMPLIANCE**  
Higher attendance in-person or remote FU visits

**IMMEDIATE VASCULAR DETECTION**

- Patient only attends a clinical visit if a problem is detected
- Efficient video-consultancy
- Simpler surgical interventions
- Shorter ICU stays

INTELVASC



WILSON SONSINI  
 US Patent number 85597361

**STEP 1: Insertable Sensor**  
The subcutaneous sensor is inserted by the clinician in an outpatient setting and the patient is discharged

**STEP 2: At-Home Scanning**  
The patient scans the implant as indicated by the clinician

**STEP 3: Clinical Review**  
The physician can access and review the Clinical data and will receive an alert if a problem is detected

## MOBILE AND CLOUD DASHBOARD

### Mobile

Easy to download and install, full walkthrough and patient guidance

### Cloud

Central backend for data processing and metric display

### USER FOCUSED EXPERIENCE

- The mobile app is aimed to be a convenient way to collect data from the reader and upload this to the cloud
- The app will guide the patient how to use and take readings as well as display important information on scheduled checkups and offer a way to communicate with clinicians directly
- The clinical backend will display important trends on the status of the intervention and repair; the clinicians can edit the frequency that the mobile app reminds the patient to take readings

## DISCOVERING WASTED RESOURCES AND COSTS PATIENT PATHWAY

Current standard	Q1	Q2	Q3	Q4	Q5	Q6	Q7	Q8
Vascular Decline	●	●	●	●	●	●	●	●
Clinical Checkups	●	●	●	●	●	●	●	●
With Interface	●	●	●	●	●	●	●	●
Vascular Decline	●	●	●	●	●	●	●	●
Clinical Checkups	●	●	●	●	●	●	●	●
At Home Scan	●	●	●	●	●	●	●	●

## REMOTE MONITORING

**REMOTE MONITORING**

- 1 MINUTE SCAN
- COMFORT OF YOUR LIVING ROOM
- PHONE REMINDERS AND ALARMS ONLY AS OPTIMALLY NECESSARY

COR	LOE	Recommendations
2a	C-LD	8. In patients with PAD who have undergone endovascular procedures without new lower extremity signs or symptoms, it is reasonable to perform ABI and arterial duplex ultrasound surveillance within the first 1 to 3 months postprocedure, then repeat at 6 and 12 months, and then annually. <sup>11,12</sup>
2b	B-NR	9. In patients with PAD who have undergone infrain-guinal, prosthetic bypass graft(s) without new lower extremity signs or symptoms, the effectiveness of ABI and arterial duplex ultrasound surveillance is uncertain. <sup>13,14,15,16</sup>
2a	C-LD	10. For patients with PAD, telehealth can be an alternative mode for vascular evaluation and management and longitudinal follow-up, but the use of these visits should be consistent with the urgency of presenting symptoms. <sup>17</sup>

**Proactive Disease Management:**

Transition from reactive approaches addressing acute complications (e.g. thrombosis) to strategies that anticipate and prevent events in vascular care.

**Enhanced Surgical Focus:**

Allow surgeons to shift from high-risk, repeat procedures to more complex cases, improving outcomes for patients with vascular disease.

**Improved Patient Compliance:**

Increased adherence to care protocols due to the convenience of at-home scanning.

**INTELVASC**<sup>TM</sup>  
VASCULAR INSIGHTS

*Experimental Data*

## PROTOTYPE DEVELOPMENT

Left: Free photoelectric ultrasonic array prototype. The array has 8 elements and is 15.3 mm in total. The central frequency is 4 MHz. Right: Comparison of flow velocity measurements performed at the same station with the prototype and with the L11-5v probe. The error bars demonstrate the standard deviation over 50 measurements.

Graph showing Flow velocity (m/s) vs Distance from occlusion (cm). Legend: Low grade stenosis (30%), High grade stenosis (80%).

Prototype: 8 elements, 15.3 mm length  
L11-5v Probe: 138 elements, 38.4 mm length

## HIGHLY RELIABLE TECHNOLOGY - FLOW MODEL RESULTS

\*Results obtained at our engineering labs at NCSU using a Venoscans Research Scanner by Feb 2024

Graph showing Flow score vs Flow rate setting (AU).  $R = 0.99$

Our ultrasonic sensor could detect flow reliably at **99% Sensitivity** relative to known inflow

### First cadaver study

Insertion of the sensor prototype over the perfused artery in a cadaveric leg. The sensor was connected and controlled by the Verasonics scanner

13

### Spatial variations of flow can be attributed to

- Turbulence or Very high flow velocity (both due to stenosis)

14

### PATHWAYS OF EXPLORATION DEVELOPMENT PATHS

**AT HOME FLOW ASSESSMENTS FOR SENSIBLE MONITORING**

**NON-INVASIVE WEARABLE**

**NON-INVASIVE IMPLANTABLE**

**NON-INVASIVE PASSIVE ASSESSMENTS**

- USER FRIENDLY
- CONVENIENT + NO MAINTENANCE
- CHARGEABLE

**WEARABLE**

**IMPLANTABLE**

**ON DEMAND, DIRECT FLOW ASSESSMENT**

- REPEATABLE AND RELIABLE MEASUREMENTS
- VERSATILITY ON MANY OTHER VESSELS
- ONE TIME IMPLANT DURING INITIAL SURGERY

15

### MEET THE TEAM

#### EXECUTIVE TEAM

- ANDREW ORTEGA-VERGADUER**  
Chief Executive Officer  
15y Commercial, Product launch Medical Devices
- KYLE FREDERICKS**  
Chief Technical Officer  
Former engineering director World's 1st wearable ultrasound
- MARIE MILLER, PhD**  
Chief Scientific Officer  
Professor Biomedical Engineering PhD Physics Ultrasound for Medical Applications
- JEAN BISHMUTH, MD**  
Chief Medical Officer  
Chief, Division of Vascular Surgery USF/FGH Medical Director of HW USF/FGH Professor of Surgery

#### ADVISORY BOARD

- ALICK LIU**  
Chief Product Officer  
20 y biomedical engineering Development Medtronic's pacemaker
- MARK ADELMAN, MD**  
Medical Chief Emeritus  
Prof. Vascular Surgery, NYU Langone Scientific Board VETH Symposium
- JAMES PRAY**  
Executive Advisor  
Former President at Penumbra MedTech Advisor & Investor
- DAVID UFFER**  
Executive Advisor  
25y experience Medical Device & Investor

- DAVID KURAJUNTILA**  
Executive Advisor  
CEO at Aio
- FRANK VETH, MD**  
Medical Advisor  
Prof. Vascular Surgery Cleveland Clinic President VETH Symposium
- HENRIK SILJESSEN, MD**  
Medical Advisor  
Professor of Vascular Surgery Rigshospitalet Denmark Former President ESVS
- JAY MATTHEW, MD**  
Medical Advisor  
Interventional Cardiology MedStar Memorial Hospital Scientific Committee VIVA Symposium

**Corporations**

**Educational Institutions**

**Medical Institutions**

**Scientific Societies**

## Monitoring Vascular Stenosis, Anytime, Anywhere.

Early Diagnosis Remotely From Home

US Patent number 63/597,301 Submitted by **WILSON SONSSINI**

**Vascular Sensor**

Monitoring vessel stenosis to anticipate vascular events

**Wirelessly Powered**

An Ila biosensor. Powered only when needed

**Easy To Implant, Easy To Use**

Surgeon focused design, patient focused experience

**AI Powered - Reliable Metrics**

A powerful data backed assisting informed medical decisions