

Improved Fiberoptic Based Imaging With The FORS Lumiguide System:

Fiberoptic Sensors Within A Guidewire & A Special 3D Hub Allow Imaging Guidance With Any Catheter: How Does It Work And Experience

Carlos H. Timaran, MD
Chief, Endovascular Surgery
Sam H. Phillips Jr MD Distinguished Chair in Surgery
Professor of Surgery

Division of Vascular and
Endovascular Surgery

Department of Surgery

UTSouthwestern
Medical Center


Disclosures

- Honoraria / Consultant / Research
 - Cook Medical Inc
 - W. L. Gore & Assoc
 - Philips Healthcare
- Some devices presented here are investigational and have not been approved by the FDA

UTSouthwestern
Medical Center

Fiber Optic RealShape (FORS)

- FORS: enables real-time 3D visualization of the full shape of devices inside the body, without the need for radiation




UTSouthwestern
Medical Center

Fiber Optic RealShape (FORS)




UTSouthwestern
Medical Center

Fiber Optic RealShape (FORS)




UTSouthwestern
Medical Center

Fiber Optic RealShape (FORS)



UTSouthwestern
Medical Center

Fiber Optic RealShape (FORS)




UTSouthwestern
Medical Center

What is New?


LumiGuide, powered by FORS technology

LumiGuide equipment



LumiGuide Docking Base LumiGuide Engine

FORS-enabled single-use devices



FORS-enabled Guidewire

3D Hub



Docking Top


UTSouthwestern
Medical Center

LumiGuide, powered by FORS technology

- Works with the catheters you're familiar with, thanks to **3D Hub**
- Has AI-based semi-automatic registration
- Visualizes tight bends such as SOS-like catheters and includes estimates for degraded sections



3D Hub





Division of Vascular and
Endovascular Surgery

UTSouthwestern
Medical Center

LumiGuide wire use in USA Pending CE mark

- In body section from 120cm to 160cm
 - Better handling of tethered device
 - Support longer (up to 135cm) catheters



UTSouthwestern
Medical Center

3D catheter agnostic guidance With LumiGuide

Connecting the 3D Hub to catheter of choice



Introducing the FORS guidewire through the 3D Hub and catheter of choice



Catheter selection and registration



3D Guidance for both FORS Guidewire and compatible Catheter



UTSouthwestern
Medical Center

First experience with LumiGuide – 3D Hub

- 82yr old female with asymptomatic 51 mm complex AAA
- Planned for CMD 4 vessel
- TV were cannulated with FORS




Division of Vascular and
Endovascular Surgery

UTSouthwestern
Medical Center

First experience with LumiGuide – 3D Hub

- 73yr old female with asymptomatic 55 mm complex AAA
- Planned for CMD 4 vessel
- TV were cannulated using LumiGuide and 3D Hub



TV cannulated with FORS using multiplanar visualization and minimum fluoroscopy

Division of Vascular and Endovascular Surgery

UT Southwestern Medical Center



LumiGuide using UE access

84-year-old male

Chief Complaint: Symptomatic 62 mm Juxtarenal aneurysm

Comorbidities: HTN, HLD and vertigo

Past medical history: No family history of aneurysm









Type III Arch





Treatment Plan

- Branched endovascular aortic aneurysm repair with the Gore® Excluder® Thoracoabdominal Branch Endoprosthesis (TAMBE)







Early experience with LumiGuide

Variable	Procedure with FORS (n=39)
FORS cannulation technical success	81.7%
Total procedure duration, min	440
Fluoroscopy time, min	82.36
Total DAP, Gy \cdot cm 2	215.98
Total RAK, mGy	1661.02

Department of Surgery
UTSouthwestern Medical Center

Conclusion

- FORS allows 3D multiplanar navigation with AI semi-automatic registration, which enhances target vessel catheterization with acceptable technical success and potentially reduced radiation
- LumiGuide allows the use of any catheter through both upper extremity and transfemoral access
- The 160-cm LumiGuide wire provides improved torqueability and navigation
- FORS is becoming a gamechanger for endovascular procedures due to the multiple benefits and applications

Division of Vascular and Endovascular Surgery
UTSouthwestern Medical Center

